

DATA SHEET

# Adapter TP-NeoRec16-ECG12

Electrode adapter cable for connecting up to 17 any electrodes with TouchProof 1.5 mm connectors to NeoRec 16 amplifier.

**REF**

**003-5-224**

**INTENDED USE**

For connecting EEG/ECG electrodes to an amplifier.

**SET**

- Adapter TP-NeoRec16-ECG12.



**DESCRIPTION**


Adapter TP-NeoRec16-ECG12 is a passive adapter for connecting EEG or ECG electrodes with individual connectors TouchProof 1.5 mm (DIN 42 802-ST) to EEG amplifier NeoRec 16.

**SPECIFICATION**

Connector 1 type (to cap)	TouchProof 1.5 mm DIN 42 802-BU x 17
Connector 2 type (to ampifier)	ST40X-18S
Maximum number of electrodes to connect	17
Electrode positions	EEG: FP1, FP2, O1, O2, F3, F4, C3, C4, P3, P4, T7, T8, Fz, Cz, Pz, GND, AUX ECG: R, L, F, N, C1, C2, C3, C4, C5, C6
Adapter / Extension lenght	0.75 m
Net weight	70
Gross weight	75

## TP-NEOREC16-ECG12 MOD.001

Pin layouts of common connector ST40X-18S.  
(Manufacturer's code)



#76-01M17