

DATA SHEET

Adapter NeoRec21-DB25F

Adapter for connecting the cap with a DB-25M common connector to an amplifier NeoRec 21.

REF

003-5-090

INTENDED USE

For connecting the electrode cap the amplifier.

SET

- Adapter NeoRec21-DB25F.



DESCRIPTION

Adapter NeoRec21-DB25F is a passive adapter for connecting an EEG cap with a DB-25M common connector to the EEG amplifier NeoRec 21.

SPECIFICATION

Connector 1 type (to cap)	DB-25F
Connector 2 type (to amplifier)	ST40X-24S
Maximum number of electrodes to connect	22
Adapter / Extension length	0.3 m